



Efficient And Versatile CRISPR/Cas9 Platform Facilitates Precise Genetic Modification In Mammalian Cell Lines


Huanyu Jin, Andreia Sommer, Jingyuan Cao, Vladimir Pak, Hiral Dantara, Yin Zhang, Zhaoti Wang, Padmaja Tummala, Nga Nguyen, Diana Nguyen, Jinling Li, Ruby Yanru Chen-Tsai

Applied StemCell Inc., Milpitas, CA 95035


Introduction

- CRISPR/Cas9 technology has revolutionized genome engineering and has emerged as a reliable and versatile tool for genome editing in various organisms and cells.
- Precise mutations, gene disruption, mutation corrections and insertions has enabled better understanding of the genetics and mechanism of diseases.
- CRISPR gene editing in cell lines, the workhorse of preclinical and biomedical research, enables the generation of unlimited *in vitro* models with precise gene modifications and advanced gene expression design that are physiologically relevant.
- To accurately model diseases, as well as for precise single or bi-allelic manipulations, there is a strong dependency on homology directed repair (HDR) and gRNA selection strategies.
- Here, we present data demonstrating the versatility of the technology in editing cell lines:
 - Editing a variety of mammalian cell lines including hard-to-transfect blood lineage cells and stem cells
 - Case studies highlighting complex modifications such as double gene knockout, large knock-in of transgenes into specific endogenous locus and safe harbor locus; gRNA selection strategies to ensure specific mono- and bi-allelic modifications.
 - Efficiency of generating knock-in, knockout and point mutations in various cell lines, including pluripotent stem cells.


Knock-in of a Large "Fusion" Transgene into a Specific Locus in a Human Cancer Cell Line



Large Transgene into hH11 Safe Harbor Locus in Induced Pluripotent Stem Cells



Cas9 Mediated Genome Engineering




Cell Lines Amenable to CRISPR/Cas9 Gene Editing*


Human		Other Species			
Cancer Cells	786-O (Renal carcinoma) A-375 (melanoma) Breast cancer cells CHLA-10 (Neuro-ectodermal) DLD-1 (Colon cancer) Gist-T1 (Gastrointestinal tumor) HKT 116 (Colon cancer) HEK293 (Kidney) HEK293T (Kidney) HEK293Brfl HeLa (Cervical cancer) HepG2 (Liver sarcoma) HT1080 (Fibrosarcoma) HT29 (Colon cancer) K562 (Myeloid leukemia) KM12-Luc (Colon carcinoma) LnCaP (Prostate cancer) MDA-MB231 (Breast cancer) NCI-H2228 (Adenocarcinoma) RKO (Colon cancer) TC32 (neuro-ectodermal cancer) SH-SY5Y (Glioblastoma)	Epithelial Cells	A-549 (lung) BEAS-2B (Bronchial) BT1474 (mammary gland) Follicular thyroid cell line HaCaT (Keratinocyte) HBE (Bronchial) Huh7 (Liver) MCF-10A (mammary gland) OCCM-30 (Cementoblasts) RPE-1 (retinal epithelium) SK-MEL-31 (skin epithelium) U-2 OS (bone)	Mouse	3617 (Adenocarcinoma) 4T1 (Mammary cancer) cTEC (Cortical thymic epithelial) C2C12 (Myoblasts) GD25 (Fibroblasts) Embryonic stem cells (mESCs; C57BL/6) Induced pluripotent cells (iPSC) OCCM (Cementoblasts) Tonsil epithelial cells
Stem Cells	iPSC (Healthy and various disease models) Human multipotent adult progenitor stem cells Neural stem cells (Adult and fetal) Limbal stem cells	Blood-derived Cells	BCiW-1 (Bone marrow) H929 (Bone marrow) Jurkat (T lymphocyte) K562 (Erythroleukemia) KHYG-1 (NK cell leukemia) MM.1s (B lymphocyte) MWCL-1 (lymphoma) NCI-H929 (Bone marrow) TF-1 (Leukemia) T2 cells (Lymphocyte)	Rat	CWSV-1 Cells (Hepatocytes) DAC8 embryonic stem cells (Dark agouti male) Immortalized Keratinocytes L77 embryonic stem cells (Fisher) PC13 (Follicular thyroid cells) Chondrosarcoma cells
Fibroblasts & Adipocytes	Immortalized fibroblasts U2 OS SGBS (preadipocyte/adipocyte)	Primate	COS-7 (Kidney cells) GL37 (Kidney cells) Macaque ES cells	Hamster	CHO-K1 (Ovarian) CHO-S (Ovarian)
		Insect & Others	Mosquito cells Wolf keratinocytes		

* Selected list of cell lines successfully modified in Applied StemCell's Cell Biology Lab; CRISPR cell line editing is not limited to these cell lines.

Generation of Double Gene Knockout in Hard-to-Engineer Jurkat Cell Lines



Targeted Heterozygous and Homozygous Point Mutations in Human Embryonic Stem Cells



Efficiency of CRISPR-Mediated Genome Editing in Cell Lines

Types of Projects	% Efficiency Homozygous	% Efficiency Heterozygous	% Total Efficiency	# of Projects
Knock-in	13.51	40.91	40.35	17-28
Knockout	14.38	Not Applicable	14.38	23
Point Mutation	8.73	6.69	10.99	19-37

The efficiency of CRISPR-mediated gene modifications was calculated as the number of heterozygous or homozygous clones normalized to the total number of clones screened. The total efficiency denotes the total number of clones (heterozygous and homozygous clones) identified versus the total number of clones screened.

Conclusions

- CRISPR/Cas9 technology is a versatile gene editing technology and can be used for modifying a variety of cell lines including hard-to-transfect blood lineage cells such as Jurkat, bone marrow cell lines, and pluripotent stem cells.
- CRISPR can also be used to efficiently and precisely modify genes: knockout/ disruption, point mutations and transgene knock-in. As well, with proper gRNA and targeting design, off-target activity can be avoided or minimized.
- CRISPR technology also enables the generation of complex genetically engineered cell line models, double knockout, and inducible gene expression models, in addition to being easily manipulated to generate specific bi-allelic and mono-allelic clones.
- In summary, CRISPR gene editing in cell lines provides an unlimited source of physiologically relevant *in vitro* models for basic research, drug target discovery and initial stage drug screening, and has tremendous potential for cell replacement therapeutic applications.