



Disease Modeling and Drug Screening in Neurological Disorders Via Novel iPSC-based Technologies


Huanyu Jin, Neeraja Parameswaran, Ying Pei, Xianmin Zeng

Applied StemCell Inc., Milpitas, CA 95035

Introduction

- Human iPSC technology offers the benefits of a cell line, coupled with the advantages of using human primary cells.
 - Human disease mutations can be captured in a stable cell population.
 - iPSCs can be terminally differentiated into multiple cell lineages and genetically engineered generating cell line models with the same allelic background.
- iPSC technology and its differentiation into neuronal lineage cells has benefited research in neuroscience and neurological disorders.
- Cell-based assays using iPSCs and differentiated cells have been approved by international drug regulatory agencies for neurotoxicity screening of drug candidates.
- To aid in furthering drug development and screening for PD, we have generated a panel of iPSC lines & terminally differentiated them into neural lineage cells for neurotoxicity assays and disease modeling applications.
- We describe the utility of these lines for neurotoxicity assays, including assays to determine the specificity of different neural cell types for a small range of chemicals and drugs from the Tox21 library, as well as for neuroprotective assays with dopaminergic neurons.

Experimental Design



Neurological Disease Modeling and Drug Screening With Three Panels of iPSCs


- Control Lines:** Well-characterized, integration-free control iPSC lines generated from male and female CD34+ cells (cord blood) using episomal vectors. These lines were also used for engineering isogenic lines with disease mutations and reporter lines.
- (A) Engineered isogenic iPSCs to model diseases:** From parental control iPSCs with knockout mutations of genes associated with neurological disorders.
- (B) Engineered reporter lines:** Isogenic lineage-specific reporter lines engineered from parental control iPSCs with knock-in of reporters under control of endogenous neuronal promoters; and reporters inserted in safe-harbor locus under control of lineage-specific or ubiquitous promoters.

Isogenic knock-out lines	Disease
PARK2 -/-	PD
PARK7 -/-	PD
PINK1 -/-	PD
LRRK2 -/-	PD
Park2-/-; Park7-/-	PD
Park2-/-; Pink1-/-	PD
APOE -/-	Alzheimer's disease
SOD1 -/-	ALS
DICS1 -/-	Schizophrenia
CNTNAP2 -/-	Autism
BDNF -/-	CNS

Knock-in neural lineage-specific reporters	Description
MAP2-Nanoluc-Halotag KI	Neuron reporter
GFAP-Nanoluc-Halotag KI	Astrocyte reporter
MBP-Nanoluc-Halotag KI	Oligodendrocyte reporter

Safe-harbor knock-in lines	Description
CAG-GFP, AAVS/Chr19	Ubiquitous reporter
DCX-GFP	Neuron reporter

Differentiation to Isogenic Panels of Neurons & Glia Using Neural Stem Cells as a Stable Intermediate



Characterization of Differentiated Astrocytes, Cortical & Dopaminergic Neurons

A

Immunocytochemical characterization of biomarkers expressed by differentiated astrocytes and neurons.

B

SYMBOL	NSC	NEURONS	ASTROCYTES	DESCRIPTION
ATP10B	2	1018	1026	
DCX	467	1918	182	
GREM2	11	181	182	
LHX1	15	1828	16	
LUC135508	284	490	317	
MAP2K12	2	869	7	
MAP6	1266	4310	1004	
MYT1	70	851	578	
NEUROD1	53	138	190	
NEUROG1	34	185	137	NEURONS
NEUROG2	207	6917	131	
PORK1	245	456	96	
PLXNA2	23	780	328	
POU4F2	5	138	69	
PS21	8	11	9	
SEMSA	212	442	282	
SLC17A6	12	252	161	
SLC17A8	15	111	76	
TUBB3	2926	456	8729	

Whole genome profiling of differentiated neurons and astrocytes was used to confirm expression of lineage specific markers.

C

Neurons and astrocytes can be co-cultured for complex neuronal models to evaluate the interactions between neurons and glia.


D

Electrophysiological activity of neurons can be measured using neurons cultured on MEA plates for up to 3 months; Raster plot shows neuronal activity across 64 electrodes on day 49.


E

Single pulse current (patch clamp) recording of dopaminergic neurons derived from iPSCs indicate neurons are excitable upon injection of current.

Isogenic Knockout Lines Recapitulate Phenotypes of Abnormal Neurons Derived from Patient Lines (Ex. Park2 -/-)



Lineage-Specific Reporter Knock-in iPSC Lines for Neurotoxicity Screening



Screening Drugs for Neuroprotection & Neurotoxicity


Table 1. Drugs that were neuroprotective in iPSC and differentiated neuronal cells, and used for human clinical trials.

Neurotransmitter/ MAO Inhibitors:	Rasagiline, selegiline, nicotine, topiramate, amantadine, zonisamide, taurine
Antioxidant/ Mitochondrial Stabilizers:	Resveratrol, N-acetyl cysteine, lipoic acid, epigallocatechin gallate, creatine
Anti-Inflammatories:	Rolipram, indomethacin, 7-nitroindazole, 3-aminobenzamide, phenanthridone

Table 2. Drug that were not neuroprotective in iPSC-based models but were neuroprotective in conventional cell lines and animal models:

Neurotransmitter/ MAO Inhibitors:	Donepezil, caffeine, theophylline, pergolide, apomorphine, riluzole, pramipexole
Antioxidant/ Mitochondrial Stabilizers:	Ascorbic acid, coenzyme Q10, uric acid, folic acid, ropinirole
Anti-inflammatories:	Minocycline, estradiol, clioquinol, plicamycin

Dopaminergic (DA) neurons derived from control iPSC lines were used to evaluate neuroprotection of compounds previously shown to be neuroprotective in rodent and cell line models (Table 1 and 2), when challenged with rotenone or MPP+. Cell viability was measured using the MTT assay. Only 18 out of the compounds (Table 1) were found to be neuroprotective in these iPSC-derived DA neurons, and these same compounds have been used in human Parkinson's disease neuroprotection clinical trials.



(A) Screening for 80 compounds in the Tox21 library showed differential toxicity in isogenic iPSCs, NSCs, Neurons, and Astrocytes, using MTT assay in 96-well plates and at two doses for each compound.

Conclusions

- We have developed a panel of lines including control, and engineered isogenic and reporter lines, which provide a unique advantage for disease modeling and drug screening.
- We have established robust methods for generating neurons and glia from iPSC using a NSC gateway concept from virtually all lines.
- Human neural cultures may better mimic human neurodegenerative disorders and therefore are a more relevant model to screen for drug efficacy and toxicity.
- We have shown that iPSC-derived neuronal and glial cells can be used for modeling neurodegenerative diseases as well as for neurotoxicological and neuroprotective drug screening.